



# PTS Software 3.5 release information



- **PTS software release 3.5 highlights**
  - ZETA Faster strategy
  - 24-2C and 30-2C test fields
  - Import from HFA 3
  - Export/Import DCM OPV from PTS
  - Kinetic examination controls
  - Quick Blind Spot detection
  - SOCT Structure & Function support
    - providing VF probability results for retina map overlay
    - flexibility in matching PTS patient to SOCT
  - Other improvements



# 3.5






# 24-2C and 30-2C test fields

- **24-2C and 30-2C**
  - Variation of popular 24-2 and 30-2 test fields
  - Added extra 10 test locations in Macula 10° area
  - Better chance to avoid missing paracentral defects without need to do the 10-2 testing



# Import from HFA3

- Import of single and multiple results from HFA3
- HFA3 results are imported numerically for future potential of numerical analysis (comparison, progression)
- Imported results are available as original reports or can be reprinted in PTS format
- Supported HFA 3 output types:
  - XML+PDF
  - DCM EPDF




- Exports PTS results to Dicom SOP Class for OPV modality: Ophthalmic Visual Field Static Perimetry Measurements Storage, UID 1.2.840.10008.5.1.4.1.1.80.1,
- Open and universal format for EMR data exchange
- Dicom OPV with PTS 920/925W/2000 results can be exported and imported to the PTS software

Tag	Attribute Name	Def VR	VM	Values
(0008,0005)	Specific Character Set	CS	1	ISO_IR 192
(0008,0012)	Instance Creation Date	DA	1	20200813
(0008,0013)	Instance Creation Time	TM	1	144303
(0008,0016)	SOP Class UID	UI	1	1.2.840.10008.5.1.4.1.1.80.1
(0008,0018)	SOP Instance UID	UI	1	1.2.276.0.7230010.3.1.4.3012323659.15828.1597322583.866
(0008,0020)	Study Date	DA	1	20191213
(0008,0021)	Series Date	DA	1	20191213
(0008,0023)	Content Date	DA	1	20191213
(0008,002a)	Acquisition Datetime	DT	1	20191213162836.296
(0008,0030)	Study Time	TM	1	162836.296
(0008,0031)	Series Time	TM	1	162836.296
(0008,0033)	Content Time	TM	1	162836.296
(0008,0050)	Accession Number	SH	0	
(0008,0060)	Modality	CS	1	OPV
(0008,0064)	Conversion Type	CS	1	WSD
(0008,0070)	Manufacturer	LO	1	Optopol Technology
(0008,0080)	Institution Name	LO	1	Optopol Technology
(0008,0081)	Institution Address	ST	1	Zabia 4242-400 ZawierciePoland
(0008,0090)	Referring Physician's Name	PN	0	
(0008,1050)	Performing Physician's Name	PN	1	perimeter...
(0008,1070)	Operators' Name	PN	1	perimeter...
(0008,1090)	Manufacturer's Model Name	LO	1	PTS 2000
(0010,0010)	Patient's Name	PN	1	1*ks
(0010,0020)	Patient ID	LO	0	
(0010,0030)	Patient's Birth Date	DA	1	19831220
(0010,0040)	Patient's Sex	CS	1	O
(0010,1010)	Patient's Age	AS	1	035Y
(0010,1040)	Patient's Address	LO	0	
(0010,4000)	Patient Comments	LT	0	
(0018,0015)	Body Part Examined	CS	1	EYE
(0018,1000)	Device Serial Number	LO	1	7700001
(0018,1020)	Software Version(s)	LO	1	SW:3.5.0.3k:FW:20180904
(0020,000d)	Study Instance UID	UI	1	1.2.276.0.7230010.3.1.2.3012323659.15828.1597322583.864
(0020,000e)	Series Instance UID	UI	1	1.2.276.0.7230010.3.1.3.3012323659.15828.1597322583.865




# Kinetic examination controls

- Manual isopters drawing and modification for improved control over the kinetic results
- Assigning custom colors to isopters
- Retesting a vector with a single click in examination interface



# Quick Blind Spot detection

- New method of finding Blind Spot with use of 1 – 27 test locations
- More flexible - better coverage of potential blind spot areas (18°x 6°) – works in cases when „classic” method fails
- Faster identification– often just 3 exposures are enough to identify the Blind Spot (vs 12 in „classic” method)
- Available only for PTS 2000





- **Providing VF numerical results with probability results for retina map overlay**
  - The pattern deviation probability symbols are exported with a location within VF and can be presented as an overlay to the retinal NFL thickness or significance
- **Added flexibility in matching PTS patient to SOCT**
  - If the name or surname are not ideally matching or are swapped, SOCT will present list of similar patient records form PTS database

